



Orthorobot Wire Equivalents

Trademark	Company	Orthorobot-Equivalent
Align	Ormco	Ni-Ti
Align TS	Ormco	Heat active Ni-Ti
Bendalloy	American Ortho	Beta Titanium
Bendaloy	Rocky Mountain	Beta Titanium
Beta – T		Beta Titanium
Beta – Titanium	Acme Monaco	Beta Titanium
Beatblue	Highland Metals	Beta Titanium
Betaflex Titan	Forestadent	Beta Titanium
Beta Force	Ortho Tech	Beta Titanium
Beta III	3M Unitek	Beta Titanium
BioStarter Martensitic	Forestadent	Heat active Ni-Ti
Black Diamond	Class One	Ni-Ti
CNA	Ortho Organizers	Beta Titanium
Copper Niti	Ormco	Ni-Ti
Equire	Dentaurum	Heat active Ni-Ti
Heat-Activated	3M Unitek	Heat active Ni-Ti
NEO-Sentalloy	GAC	Heat active Ni-Ti
Nickel – Ti		Ni-Ti
Noibium	Ormco	Beta Titanium
Nitanium	Ortho Organizers	Ni-Ti
Reflex	TP Ortho	Ni-Ti
Rematitan Lite	Dentaurum	Ni-Ti
Rematitan Special	Dentaurum	Beta Titanium
Resolve	GAC	Beta Titanium
Sentalloy	GAC	Ni-Ti
Tensic	Dentaurum	Heat active Ni-Ti
Timolium	TP Ortho	Beta Titanium
Titanal	Lancer	Ni-Ti
Titanol SE, Triple F, Low F	Forestadent	Ni-Ti
Titanol Martensitic	Forestadent	Heat active Ni-Ti
Titanium Molybdenum	Acme Monaco	Beta Titanium
TMA	Ormco	Beta Titanium
Tru Flex	Ortho Tech	Ni-Ti
Vision SE	Nexadental	Ni-Ti
Vision Thermal	Nexadental	Heat active Ni-Ti

All trademarks and registered trademarks are property of their respective owners.